



Healthcare



# CUVIS

JOINT



**UP**

- Safety
- Consistency
- Surgical Accuracy
- Flexible



**DOWN**

- Surgical Error
- Complications
- Side Effect
- Possibility of Re-operation

# Overview

CUVIS Joint® robot system for artificial joint surgery is the most advanced surgical equipment capable of 3D pre-planning, virtual surgery and precise cutting to provide accurate and precise surgery results.

It is a user-friendly surgical robot equipped with the values of Innovative, Flexibility and Easiness of use with the basic offering of accuracy and safety at its best.

## 01 Main Console



- ✓ System Monitor
- ✓ Optical tracking system (OTS)

## 02 Robotic Arm



- ✓ Operating Software
- ✓ Main Controller
- ✓ Robot Arm
- ✓ Milling Tool
- ✓ Irrigation

## 03 Planner



- ✓ Surgical Planning Software

## Key Features



### Product & Class

Automation Robot surgical system (Class 3)

### Planning

3D CT Base surgical planning

### Cutting

Fully automatic cutting (Milling)

No cutting guides (jigs) required

Tibial cut including keel preparation

Fully finished femur with all the cuts and peg holes

### Safety

Emergency stop & Manual guide

### Detection

OTS (Optical Tracking System)

### Robot cutting

Max. 50mm/sec

### Robot precision

Repeat precision < 0.5mm / Positioning accuracy < 1mm

### Other features

Wide surgical space and high freedom (6axis articulated robot)

Patient-Robot Position Guide (Surgical space check)

Speed control during cutting

Compact hardware, minimize product space

## Why Robotic Artificial Joint Surgery is good?

### 01 Personalized pre-planning

As each person has a different face, the shape of bone is also different. Meril's artificial joint surgical robot CUVIS Joint<sup>®</sup> shows the patient's bone in 3D images, and the doctor can use those images for pre-planning of surgery personalized for the patient.



### 02 Pre-selection of artificial joint and precise insertion of artificial joint

What's as important as the precise surgical plan is to select and insert the personalized artificial joint. The doctor uses robot to select an artificial joint for the patient and insert it accurately.



## Why Robotic Artificial Joint Surgery is good?

03

### Precise cutting for sub-millimeter accuracy and optimum alignment


Precise cutting serves the optimum result. CUVIS Joint® provides the correct alignment of a patient's leg axis with the sub-millimeter dimensional accuracy and precise cutting for the optimal surgical outcome.



04

### Reduction of side effect and reoperation

CUVIS Joint® reduces side effects like inequality of limb length, pulmonary embolism, and fracture. The risk of infection is also reduced because of fewer instruments in use than in conventional surgery.



# Cuvis Joint Knee Replacement Process



## STEP 1.


Patients can decide upon robotic artificial joint surgery after consulting with the doctor

## STEP 2.

CT scanning.

## STEP 3.

The scanned CT image is converted into a 3D image for diagnosing the patient's condition and make a surgical plan as required.



\*FMA : Femoral Mechanical Axis/TMA : Tibia Mechanical Axis

## STEP 4.

The patient is connected to the robot and stabilized for surgery. In the next stage, the doctor performs registration process to verify if the 3D image of the patient matches the original surgery site.

After registration process, robot reviews the data and cuts the bone precisely with respect to size, position, angle and direction of the implant decided during pre-surgery planning stage.

## STEP 5.

Insert and fix the decided implant for surgery conclusion.

## Advantages of Cuvvis Joint



### Simplicity

- › Bone registration using probe
- › Improved workspace check
- › 3D bone model generation with fast and easy CT image data
- › UI design considering user convenience
- › Reduce the surgery preparation time- 3D modeling, Non-Sterile/Sterile Diagnosis



### Flexibility

- › Various cutting options
  - Full and partial cutting
  - change cutting order
- › Intra-operative Gap Check
  - Pre/Intra/Post
- › Plan changing, Gap Balancing possible



### Safety

- › Real-time system monitoring  
Emergency stop & force freeze  
Manual guide of robot arm




### Accuracy

- › Precise pre-surgical planning executed every time
- › Sub-millimeter dimensional accuracy
- › Optimal Alignment




# Surgical Process

## 1. Scan




- CT Scan

## 2. Surgical Planning




- CT Data Loading  
- Set up surgical Plan


## 3. System Diagnosis



- Non-Sterile / Sterile  
- Tool and Base Marker Positioning




## 4. Registration




- Patient setup &  
-Registration with Patient Bone Model

## 6. Surgery Result / Implant Insertion



- Check cutting results &  
Implant insertion  
- Robot out


## 5. Cutting / Gap Check



- Measured Resection / Modified gap technique  
- Pre / Intra / Post - Resection gap check

## Intraoperative Assistance


### Intra-Resection Gap Check & Plan Changing




## References

1. Jacofsky DJ, Allen M, Robotics in Arthroplasty: A Comprehensive Review, J Arthroplasty, 2016 Oct;31(10):2353-63
2. Philip Pastides, Dinesh Nathwani, The role of newer technologies in knee arthroplasty, Orthopaedics and Trauma, 2017 Feb;31:47-52
3. Eun Kyoo Song, Jong Keun Seon, Computer Assisted Orthopedic Surgery in TKA, Recent Advances in Hip and Knee Arthroplasty, 2012:421-442
4. B. Kayani, F. S. Haddad, Robotic total knee arthroplasty – clinical outcomes and directions for future research, 2019 Oct;8(10):438-442
5. Nikhil Agarwal, Kendrick To, Stephen McDonnell, Wasim Khan, Clinical and Radiological Outcomes in Robotic-Assisted Total Knee Arthroplasty: A Systematic Review and Meta-Analysis, 2020 March; 1-17
6. B. kayani, S. Konan, J. Tahmassebi, J. R. T. Pietrzak, F. S. Haddad, Robotic-arm assisted total knee arthroplasty is associated with improved early functional recovery and reduced time to hospital discharge compared with conventional jig-based total knee arthroplasty, 2018 July, 100-B:930–7

### Site Design



BIRD VIEW



TOP VIEW

www.merillife.com

**CUVIS**  
JOINT

**Corporate Headquarters, Vapi**

Meril Healthcare Pvt. Ltd.  
H1-H3, Meril Park  
Survey No. 135/2/B & 174/2  
Muktanand Marg, Chala  
Vapi 396 191. Gujarat. India.  
T: +91 260 3052 100  
F: +91 260 3052 125  
E: orthopedics@merillife.com

**International Sales & Marketing**

Meril Healthcare Pvt. Ltd.  
612, B-Wing, Bonanaza  
Sahar Plaza, Andheri East  
Mumbai 400 059  
Maharashtra. India.  
T: +91 22 3935 0700  
F: +91 22 4047 9717

**India Sales & Marketing**

Meril Healthcare Pvt. Ltd.  
512, Midas, Sahar Plaza, J.B. Nagar  
Andheri East, Mumbai 400 059  
Maharashtra. India.  
T: +91 22 4047 9797

**Meril South America**

DOC MED LTDA  
1079 – Cep: 04077-003 - Moema  
Sao Paulo, Brazil  
T/F: +55 11 3624 5935/6

**Meril Germany**

Bornheimer Strasse 135 – 137  
D – 53119 Bonn, Germany  
T/F: +49 228 7100 4000/1

**Meril Turkey**

Imalat Ve Ticaret A.S.  
6, Mimar Sinan Mah.,  
Cavusbasi Cad, Ozde Sok  
Aydin Eksi Is Merkezi Kat:1  
Cekmekoy/Istanbul, Turkey  
T: +90 53 2272 5172



MH/BRC/KS/005 Rev. 00

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Meril does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery. The information presented is intended to demonstrate the breadth of Meril's product offerings.

A surgeon must always refer to the package insert, product label and/or instructions for use before using any of Meril's products. Please contact your sales representative if you have questions about the availability of any of Meril's products in your area.