


Evermine 50 EES-1 study

Study Highlights

- Principal Investigator: Dr. K. Sivaprasad
- Evermine 50™ EES-1 is a prospective, post-marketing, multi-center, study to evaluate safety and performance of ultrathin strut biodegradable polymer-coated with everolimus-eluting coronary stent (EES) in patients with coronary artery disease
- Safety and performance endpoints of the stent will be evaluated till 2 years
- The study demonstrated favourable safety and performance of ultra-thin biodegradable polymer EES with ultrathin strut thickness for treatment of coronary artery disease



❖ Study Design

A prospective, multicentre, post-marketing study



A total of 114 patients were enrolled




9 investigational sites across India



Clinical follow-up for 1-month and 9-months and Telephonic/clinical follow-up for 6-month, 12-month and 24-month



Angiographic follow-up at 9 months
Analysed by CBCC Global Research LLP, Ahmedabad, India



OCT follow-up at 6 months
Analysed by CBCC Global Research LLP, Ahmedabad, India

❖ Study Results




Figure 1: Summary of diseased vessel in patients treated with Evermine 50 EES, (%)




Figure 2: Cardiac status of patients treated with Evermine 50 EES, (%)

❖ QCA Analysis




Figure 3: Late lumen loss at 9-month follow-up

❖ References

1. Clinical Trial number [CTRI/2017/02/007781]

<http://ctri.nic.in/Clinicaltrials/pmaindet2.php?trialid=16395&EncHid=&userName=CTRI/2017/02/007781>

2. K. Sivaprasad et al. Initial Experience With Ultrathin Strut Biodegradable Polymer-Coated Everolimus-Eluting Coronary Stent in Coronary Artery Disease Patients. *Journal of the American College of Cardiology*. 2018, 72(13). Supplement/B220.2.

3. K. Sivaprasad, One-Year Clinical Outcomes of World's Thinnest (50 μ m) Biodegradable Polymer-Coated Everolimus-Eluting Coronary Stent in Patients with Coronary Artery Disease. TCT 2019. (25-30 September)

4. K. Sivaprasad, Clinical outcomes of world's thinnest (50 μ m) biodegradable polymer-coated everolimus-eluting coronary stent system in patients with coronary artery disease. EuroPCR 2020 (19-22 May).