

Fiona-1 Study

A clinical study to evaluate safety and performance of Fiona Levonorgestrel – Releasing Intrauterine System

Study Design

- An open label, post-marketing clinical study
- 309 subjects at 21 sites across India

Reference No.	CTRI/2016/03/006781
Study Objective	<ol style="list-style-type: none">1. The primary objective was to evaluate safety and performance of Fiona Levonorgestrel-Releasing Intrauterine system as indicated by:<ol style="list-style-type: none">a. No clinical features of intolerance or complications developing within 2 months of study system insertionb. Non-failure of contraception in presence of Fiona Levonorgestrel-Releasing Intrauterine systemc. No chronic complication up to 5 yearsd. No auto-expulsion or dislocation of the study device2. Pharmacokinetic evaluation of plasma Levonorgestrel concentration up to 5 years
Efficacy Endpoint	Proportion of failure of contraception as indicated by events of pregnancies in presence of study device or due to auto-expulsion of study device <ul style="list-style-type: none">• Number of Intra-uterine pregnancies• Number of ectopic pregnancies
Safety Endpoints	<ol style="list-style-type: none">1. Proportion of population presenting intraprocedural perforation of uterus2. Proportion of population presenting with one or more of the following conditions within 1 month following insertion of study device<ol style="list-style-type: none">a. Post procedural perforation of uterusb. Significant endometrium injuryc. Myometrium injuryd. Pelvic inflammatory disease, severe Leucorrhoea or endometritis not relieved by medical intervention and requires removal of study device.e. Expulsion or partial expulsion of the study device3. Proportion of subjects presenting with events of dislocation or expulsion of the study device, perforation of uterus, endometrial injury or ulcer, myometrial injury, ectopic pregnancy, salpingitis or similar conditions, suspected to be related with the study device, after one month up to 5 year from insertion of the study device
PK Endpoint	Demonstration of Levonorgestrel levels in blood up to 5 years.
Clinical Sites	21 sites across India
Sample Size	Total 309 subjects enrolled

Follow-Up	Clinical follow-up at 1 month, 3 months, 6 months, 1 year, 2 years, 3 years, 4 years and 5 years
Study Duration	Study starts date April 2016. Estimated completion date August 2022.

❖ **Reference:**

1. Clinical Trial Registry- India (CTRI): CTRI/2016/03/006781

<http://ctri.nic.in/Clinicaltrials/showallp.php?mid1=14673&EncHid=&userName=CTRI/2016/03/006781>