

Meril Extraction Free

Dry Swab Kit

Meril Extraction Free Dry Swab Kit is used to isolate the nucleic acid (RNA) from dry swab samples for clinical *In-vitro* diagnostics applications.



**No
VTM & RNA
Extraction
Required**

**Time
Saving**

**Cost
Effective**


Meril Extraction Free

Dry Swab Kit

Meril

Diagnostics

Procedure



Features

- Compatible with nasal swab
- Sterile swabs provided
- Easy Nucleic acid elution process
- Highly safe sample transportation
- High yield of RNA
- Stable & Reliable RNA elute, Compatible PCR, RT-PCR, qPCR, qRT-PCR and other downstream molecular biological experiments

Kit components

- Dry Swab (DS) Buffer
- Sample collection tubes
- Sterile nasal / oral swabs
- Tube stand
- Pack Insert / IFU

Stability

- Storage of DS buffer should be 2 to 8 °C
- Sample should be store at 4°C upto 48 hrs.

Product Code - MBTDSK-01

Description - Meril Extraction Free Dry Swab kit

Pack Size - 100 Tests

Meril Diagnostics Pvt. Ltd.

D1 - D3, Second Floor, Meril Park, Survey No. 135/2/B & 174/2, Muktanand Marg, Chala, Vapi - 396191, India

601, Midas, Sahar Plaza Complex, J.B Nagar, Andheri East, Mumbai - 400059

Tel: +91 22 30732400 | Fax: +91 22 39350777 | Customer Care: 1-800-266-3745 (Toll Free) | Email: diagnostics@merillife.com