



Meril
Cardiac
Solutions

Meril

MYOSISTTM

Intra Aortic Balloon Catheter Kit (8F)

Performance
You Can Count On

Technical Specifications

Balloon Volume (cc)	30	35	40
Guideline Patient Height (cm)	<152	152-162	162-183
Balloon Diameter (mm)	14.5	15	15
Balloon Length (mm)	225	227	259
Balloon Catheter Length (mm)		855	
Balloon Catheter Usable Length (mm)		710	
Catheter Profile (french)		8	
Inner Lumen ID (inch)		0.028	
Guidewire Size (inch)		0.025	
Console Compatibility	Datascope®: System 97, 98, 98XT, CS100, CS300, CARDIOSAVE Arrow®: ACAT, AutoCAT, AutoCAT 2, AutoCAT 2 Wave, KAAT II		

Each Kit Contains

- Catheter Tray:**
1. Balloon Catheter
 2. 60cc Syringe
 3. One-way Valve

- Accessory Tray:**
1. PTFE Coated Guidewire 0.025" x 150 cm
 2. Uncoated Guidewire 0.035" x 50 cm
 3. 8F Introducer Sheath
 4. Dilator
 5. Double Tapered Dilator
 6. 18G Needle
 7. Pressure Monitoring Tubing - 2 pcs
 8. Three-way Stopcock
 9. Datascope® Adaptor Set
 10. Arrow® Adaptor

Ordering Information

Product Code	Description
MYS0830	Myosist IAB Catheter Kit 8F 30CC
MYS0835	Myosist IAB Catheter Kit 8F 35CC
MYS0840	Myosist IAB Catheter Kit 8F 40CC



Manufactured by:

Meril Life Sciences Pvt. Ltd.

Muktanand Marg, Chala, Vapi - 396 191, Gujarat, India

T : +91 260 3052 100 F : +91 260 3052 125 W : www.merillife.com

M.L. No. G/28/1273