



MIRUS™ Linear Cutter




Cuts and separates tissue between two double staggered rows



New knife with every reload




Every reload has safety lock to prevent firing over empty reload



Dual side firing with unique safety lock to prevent inadvertent pre-firing during preparation & positioning

Push button opening makes reload removal easy and quick

High strength and high quality material featuring anti-slip properties



Central alignment of two halves with rear hinge gives better ease of handling

Stapler Code	Reload Code	Number of Firings	Number of Staples	Staple Line Length	Cut Length
MLC60	MLCR60-3.8 MLCR60-4.8	8	64	64 mm	60 mm
MLC80	MLCR80-3.8 MLCR80-4.8	8	84	84 mm	80 mm
MLC100	MLCR100-3.8 MLCR100-4.8	8	104	104 mm	100 mm