

Meril

Endo-Surgery



PROFOUND™ N

Titanium

Non-Absorbable Mesh Fixation Device for Laparoscopic Hernia Repair



Confident Fixation,
Delivered

ri
More to Life


M-HRS
Meril Hernia Repair Solutions

Confident Fixation, **Delivered** with 4 benefits for superior surgical outcomes

1. Security, **Delivered**

Trusted helical design spike for strong hold

- Unsurpassed fixation strength
- Fixates in Coopers ligament and underlying structures



2. Comfort, **Delivered**

- Comfortable handle with easy to deploy trigger
- Significantly lesser deliver force is required
- Each spike deploys with equivalent force



Ergonomic Trigger

- Easy to use trigger with an audible CLICK ACOUSTIC feedback indicating that the spike is fully deployed
- Fits to a wide range of surgeons hand sizes

Confident Fixation, **Delivered** with 4 benefits for superior surgical outcomes




Superior Titanium alloy (Ti6Al4V ELI), for Superior Fixation Strength

- Excellent biocompatibility especially when in direct contact with tissues
- Superior ductility & fracture toughness for assured fixation in tough tissues & ligaments

3. Safety, **Delivered**

Needle-less fixation




- Eliminates likelihood of inadvertent needle sticks

- **20micron extra shaft thickness** –
Eliminates the risk of shaft bending especially during fixation of rectangular mesh edges

4. More Value, **Delivered**

Available in 15 & 30 spikes count providing significant savings



BROAD RANGE OF MESH FIXATION DEVICES

Meril Endo-Surgery offers a complete, compelling and efficient portfolio of fixation devices that delivers procedure specific solutions for hernia repair.

Meril Hernia Repair Solutions offers a comprehensive and innovative range of hernia solutions that presents greatest efficacy and value to patients.



PROFOUND™ A ABSORBABLE FIXATION DEVICE RANGE

Order Code	Description
MFD15A	PLGA ABSORBABLE 15 SPIKES FOR LAPAROSCOPIC HERNIA REPAIR
MFD30A	PLGA ABSORBABLE 30 SPIKES FOR LAPAROSCOPIC HERNIA REPAIR

PROFOUND™ N NON-ABSORBABLE FIXATION DEVICE RANGE

Order Code	Description
MFD15N	TITANIUM NON-ABSORBABLE 15 SPIKES FOR LAPAROSCOPIC HERNIA REPAIR
MFD30N	TITANIUM NON-ABSORBABLE 30 SPIKES FOR LAPAROSCOPIC HERNIA REPAIR

Please refer to the package insert for complete instructions, contraindications, warnings and precautions

Meril Endo-Surgery Pvt. Ltd.

601, Midas, Sahar Plaza, J.B. Nagar, Andheri (East),
Mumbai - 400 059. Maharashtra, India.
Toll Free No.: 1800 419 4433 Tel: +91 22 3073 2400
www.merillife.com

E1-E3, Meril Park, Survey No. 135/2/B & 174/2,
Muktanand Marg, Chala, Vapi- 396 191, Gujarat, India
Tel.: +91 260 3052100.
Email: enquiry.endosurgery@merillife.com